LCMS Troubleshooting Tips

Changes in Sensitivity

m/z



- If suppressing mobile phase components were used in previous method, clean MS source and flush

Undesired Fragmentation

Causes	Solutions
lon source setting too harsh	Check source temperatures are appropriate for analyte.Check ionisation voltage is appropriate for analyte.
Collision energy too high / low	- Check and optimise collision cell gas pressure and collision energy.
lon optics	- Ensure correct voltage is applied to desolvation line

Basic Steps

and QArray.

Follow these steps to isolate where the problems is.

Check the obvious explanations first and change only one thing at a time!



No Peaks

Solution Causes

LC

MS setting issue - Check the method uses appropriate MS settings for the compounds of interest. - Ensure probe is at correct distance from orifice.

- Check protrusion of capillary from probe.
- Check spray from capillary.
- Check ion source temperatures / gas flows are as
- expected and stable.
- Check collision gas pressure is correct and stable. - Ensure analyte elutes within event window set in
- method.
- Check for adducts such as Na⁺, NH_a⁺, K⁺. Consider different glassware, make up new mobile phases.
- Check for matrix effects such as ion suppression. Consider different LC methods and sample preparation procedures.
- Concentration injected is below limit of detection. - Perform and check MS tune.
- Check LC outlet tubing is connected to ion source. - No mobile phase flow, possibly purge valve left
 - open - Purge the system including injector to remove possible air bubbles in pump.
 - Purge LC system using isopropanol to ensure check valves are working correctly.
 - Check for crimped or damaged tubing.

Compounds	- Check mobile phase composition is correct.
not retained	- Check correct analytical column type is being used.
or retained	- Increase run time.
time in method	- Increase solvent strength.
conditions	- Check correct flow rate is being achieved.

LC system.

Low Pressure

Causes	Solutions
Partial leak in system	 Check all connections and retighten any which have leaks.
Flow rate	 Check the method has the correct flow rate. Test the flow rate accuracy using a calibrated flow rate meter or collect a specific volume and monitor the time required. Replace worn out or damaged pump seals.
Method	Check if method is using correct temperature and correct solvents.If a column section valve is used, check correct column selected.
Incorrect column	 Use correct column with correct dimensions and particle geometry.
Column temperature too high	 Set adequate column temperature and check no column damage if exceeded column temperature limit.
Airlock in LC tubing	 Remove tubing from degasser and ensure flow under gravity. Reconnect and purge pumps in isopropanol.
Stuck check valve	 Purge LC system using isopropanol and ensure check valves are working correctly. Sonicate the check valves in isopropanol.
Sensor malfunction	- Repair or replace pressure sensor.

Poor Mass Accuracy (HRAM Instruments)

Solutions Causes

 $\Delta m > expected$

MS out of tune - Perform and check system tune.

TOF Calibration - Perform TOF calibration.

Calibration - Ensure sample analytes are within calibration range and performed adjust if required. incorrectly

Detector - Dilute sample or adjust injection volume. saturation

Check the Basics:

- Power & electrical connections
- Communication cables
- Sample preparation
- Analytical conditions
- Mobile phase preparation
- Needle rinse & seal washes • Solvent flowing / no air bubbles
- (oil level & gas ballast) MS vacuum

Roughing pump

• Argon gas cylinder (level and pressure)

• LC pump pressures

• Ion source maintenance

• Gas generator (pressure readings)

Identify the Cause:

- Clearly define the problem
- Review sample and maintenance logs to identify trends in the data or possible problem indicators
- Use a logical sequence of events to isolate possible causes

Document Everything:

- Document typical values such as analyte retention times and normal initial LC operating pressures. Use as a benchmark to indicate deterioration in system performance
- Document all troubleshooting steps and results. This may help you identify and solve the next problem faster
- Always inject a known sample and compare to previous data as a reference to ensure restored performance

Still having problems with your instrument? Let us know at

office@shimadzu.eu.com

Carry Over	
Causes	Solutions
Inappropriate wash settings	- Check wash solution and wash settings.
Sample concentration	- Use lower concentrated sample or inject less.
Column	- Elush column: replace quard columns and analytical column

Sample issues Prepare fresh samples. Review injection volume in sequence / method.

- Ensure the sample is in the correct position in the autosampler.
- Check for sample adsorption issue.
- Check for air pockets trapped in bottom of

Sample flowing - Check divert valve settings if applicable. to waste

vial or well.



Changing Retention Times

Causes	Solutions
Flow rate	- Check the method uses the correct flow rate.
Insufficient equilibration	 The column should be equilibrated using at least 10 column volumes (check with manufacturer's manual). Allow more time or column volumes to equilibrate column between injections.
Poor temperature control	- Ensure column oven temperature is accurate.
Change in column dimension	- Ensure the correct column including dimensions are being used.
Change in column stationary phase environment	 Do not use a column which has ion pairing reagent for other mobile phases due to memory effects. Stationary phase 'de-wetted'.
Improper mobile phase	Ensure the mobile phase is accurately prepared and the correct lines are being chosen on the method.Ensure the pump is accurately dispensing mobile phase.
Instrument leaks	- Check for loose fittings throughout the system.
Air bubble in pump	- Purge pump via purge valve.
Sample diluent	- Inappropriate sample diluent for column.
Needle rinse	- Needle rinse solvent reaching the column.
Faulty or 'sticky' check valve	- Purge LC system using isopropanol to ensure check valves are working correctly.

High Pressure	
Causes	Solutions
Flow rate set too high	- Reduce flow rate setting.
Blocked column	- Backflush column (if permitted) or replace column.

Changes in MS Resolution

Solutions Causes - Perform and check tune. MS out of tune

contamination	if applicable.
Injector issue	 Changes in dispensing volume of injector, use a system suitability sample to determine volume changes. Check batch / method details to ensure the correct volume was programmed. Increase needle and loop flushing protocols to ensure no carryover from injection. Purging the injector metering pump.
LC Gradient	- Insufficient time at strong solvent conditions during gradient program. Increase based on column dimensions.

Incompatible mobile phase (precipitated buffer or immiscible)	 Use correct mobile phase; wash column and re-equilibrate.
Improper column	 Use correct column with correct dimensions and particle geometry.
Injector blockage	 Clear blockage (review needle, loop, valve assembly and HPV outlet).
Guard column / cartridge blockage	- Replace or remove guard column.
Column in-line filter blockage	- Replace or remove in-line filter.
Column temperature too low	- Set adequate column temperature.
Sensor malfunction	- Repair or replace pressure sensor.
Pump in-line filter blockage	- Replace in-line filter.
Blocked tubing	- Replace blocked tubing as necessary.



